

# **Vidar**

## **Digitizer Tools Version 2.0**

**Digitizer Diagnostics**  
**Firmware Download**  
**Passcode Editor**

---

## VIDAR Part Number 18210-001 Rev C

© Copyright 2009 by VIDAR Systems Corporation. All rights reserved.

No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior written permission of VIDAR Systems Corporation.

VIDAR will provide support as per the terms and conditions set forth in the license agreement. Other terms and conditions may have been presented by VIDAR resellers. If these terms differ from those offered by VIDAR, it is the responsibility of the reseller, not VIDAR to fulfill these obligations. Further, support will only be given in the following situations:

- Installation and modifications are made by qualified personnel,
- The equipment is used in accordance with the instructions in this manual.
- No other software or hardware other than those sanctioned by VIDAR or distributed by VIDAR with Clinical Express are used on the Clinical Express workstation.

# Contents

- About this manual..... 4
- Digitizer Diagnostics..... 5
  - Installing Digitizer Diagnostics ..... 5
  - Using Digitizer Diagnostics ..... 5
- Firmware Download ..... 10
  - Installing Firmware Download..... 10
  - Using Firmware Download..... 10
- Passcode Editor ..... 13
  - Installing Passcode Editor ..... 13
  - Using Passcode Editor ..... 13
- Troubleshooting ..... 17
- Appendix..... 23

# About this manual

---

This manual provides step-by-step instructions for installing and using these setup programs for Vidar film digitizers using **Windows XP, and Windows 7 (32 and 64 Bit)**:

- **Digitizer Diagnostics** comprehensively tests an attached digitizer.
- **Firmware Download** displays an attached digitizer's current firmware versions and assists in updating firmware.
- **Passcode Editor** enables a new passcode to be assigned to a digitizer.

# Digitizer Diagnostics

---

Digitizer Diagnostics performs comprehensive testing of an attached film digitizer. When testing is complete, you can send the results to Vidar Technical Support for analysis.


## Installing Digitizer Diagnostics

1. Download *DigitizerDiagnostic.exe* from the Vidar Web site. [www.filmdigitizer.com](http://www.filmdigitizer.com). Browse to the Support page and then browse to the Troubleshooting page. Save the file in a convenient location on the local computer.
2. Double-click *DigitizerDiagnostic.exe* to launch the installation. This process can take up to 30 seconds before you see any indication that the installation has started.
3. Follow the instructions provided by the installation shield. It may be necessary to install other prerequisite programs as part of the installation. These will automatically be installed as part of the installation process.
4. Follow the on screen instructions to complete the installation.

## Using Digitizer Diagnostics

**Note:** The film digitizer must be powered on for at least 30 minutes before running diagnostic tests.

1. From the Desktop double click > **Digitizer Diagnostics**.
2. You will see the warning below. Click **OK** to continue.




3. Wait while the application loads.




**Note:** The program will notify you if it cannot communicate with the digitizer. In this case, go to the "Troubleshooting" chapter.

5. Click the **Start** button in the lower right corner to start the diagnostics.



6. Wait while the diagnostic tests are performed.



7. When the tests are complete, click the **Contact Form** button in the lower right corner.

8. Fill in the fields in the **Contact** form, then click **Next >>**.


The screenshot shows the 'Digitizer Diagnostics 2.0' window. On the left, a sidebar lists three steps: 'System check' (checked), 'Requirements check' (checked), and 'Run digitizer tasks' (checked), with a 'Done' button below. The main area displays the 'Contact Information' form. At the top, it shows 'Model: DiagnosticPRO Advantage (USA model)' and 'Serial number: 351727'. The form title is 'Contact Information'. Below the title, it says 'Please fill out the form below and click the Next button below. Fields in bold are mandatory.' The form fields are: 'Digitizer Serial No.' (351727), 'Contact Name' (John Customer), 'Company Name' (empty), 'Address' (empty), 'Phone' (703-471-7070) with a note '(Include country and area codes)', 'Email Address' (john.customer@hospital.com), and 'Comments' (empty text area). At the bottom right, there are '<< BACK' and 'NEXT >>' buttons. The Vidar Systems Corporation logo is in the bottom left corner.

Digitizer Diagnostics generates a report of its findings. The report is located in the *My Documents\Film Digitizer\Report* folder. (A typical report is reproduced in the Appendix.)


9. Send the diagnostic report to Vidar using one of the following methods:
- If Digitizer Diagnostics detects e-mail capability on the computer, a new message window automatically appears with the diagnostic report attached. Send the message.  
or
  - If Digitizer Diagnostics does not detect e-mail capability on the computer, you can copy the report file to a computer that has e-mail capability, then send the report as an attachment to [medtech@vidar.com](mailto:medtech@vidar.com).  
or
  - Print the report and fax it to 1.703.471.7665.



10. *To reset the digitizer:* Click **OK** in the dialog shown below. Digitizer Diagnostics will automatically close when the reset is complete. *or*  
*To exit the application without resetting the digitizer:* Click **Cancel**.



**Note:** You must reset the film digitizer before attempting to digitize films.



# Firmware Download

Firmware Download determines a film digitizers current firmware versions, and enables you to update firmware if newer versions are available.


**NOTE: This should never run unless instructed to do so by a VIDAR authorized agent.**

## Installing Firmware Download

1. Download *FirmwareDownload.exe* from the Vidar FTP or PartnerNet sites. Save the file in a convenient location on the local computer.
2. Double-click *FirmwareDownload.exe* to launch the installation. This process can take up to 30 seconds before you see any indication that the installation has started.
3. Follow the instructions provided by the installation shield. It may be necessary to install other prerequisite programs as part of the installation. These will automatically be installed as part of the installation process.
4. Follow the on screen instructions to complete the installation.

## Using Firmware Download

1. From the Desktop double click on > **Firmware Download**.
2. You will see the warning below. Click **OK** to continue.




3. Wait while the application loads.



**Note:** The program will notify you if it cannot communicate with the digitizer. In this case, go to the “Troubleshooting” chapter.


4. Firmware Download first reads and displays the digitizer’s firmware versions.

- If the **Status** column shows “valid” for all firmware versions, then no further action is required. Exit Firmware Download.
- If firmware updates recommended:
  - a. Open the drop-down list and select the desired firmware package to install.
  - b. Click the **Start** button in the lower right corner.



5. Wait while the firmware is updated.

**CAUTION:** Do not turn off the digitizer or computer during the firmware updating process.



6. When firmware updating is complete, the digitizer will be automatically reset. When resetting is complete, Firmware Download will automatically close.

# Passcode Editor

Passcode Editor enables you to assign a new passcode to a film digitizer.


**NOTE:** This should never run unless instructed to do so by a VIDAR authorized agent.

## Installing Passcode Editor

1. Download *PasscodeEditor.exe* from the Vidar FTP or PartnerNet sites. Save the file in a convenient location on the local computer.
2. Double-click *PasscodeEditor.exe* to launch the installation. This process can take up to 30 seconds before you see any indication that the installation has started.
3. Follow the instructions provided by the installation shield. It may be necessary to install other prerequisite programs as part of the installation. These will automatically be installed as part of the installation process.
4. Follow the onscreen instructions to complete the installation

## Using Passcode Editor

1. **Start > All Programs > Vidar > Passcode Editor.**
2. You will see the warning below. Click **OK** to continue.




3. Wait while the application loads.



**Note:** The program will notify you if it cannot communicate with the digitizer. In this case, go to the “Troubleshooting” chapter.

4. In the **New password** field, enter the new password.



The word “invalid” will appear if you enter an unacceptable passcode.

Old password	<input type="text" value="MKLIZAMBVLFLBN"/>	
New password	<input type="text" value="RAB901BNTW0291"/>	invalid ←
Confirm new password	<input type="text"/>	

5. Enter the new passcode again in the **Confirm new password** field. If the contents of the **Confirm new password** field match the contents of the **New password** field, the word “matched” will appear.

Passcode Editor 2.0

Model: DiagnosticPRO Advantage (USA model)  
Serial number: 351727

System check  
Requirements check  
Run digitizer tasks  
Done

Click the UPDATE button below to continue.

**Warning:** Loading an incorrect passcode into the digitizer will damage the product with repair needed at VIDAR's headquarters. Only use this program if instructed to by VIDAR or your authorized reseller. This application is NOT capable of checking validity of the entire passcode. Double check the entry to avoid service issues. After the passcode is changed, the software will automatically reset the film digitizer. This is necessary and may take several minutes.

Old passcode	<input type="text" value="NGZOPBJOVFLNN"/>	
New passcode	<input type="text" value="NGZOPBJOVFLNN"/>	valid
Confirm new passcode	<input type="text" value="NGZOPBJOVFLNN"/>	matched


VIDAR systems corporation  
The Image of Reliability

UPDATE


The word “unmatched” will appear if the contents of the **Confirm new password** field do not exactly match the contents of the **New password** field. If this happens, correct the passcode confirmation.

Old password	<input type="text" value="MKLIZAMBVLFLBN"/>	
New password	<input type="text" value="MKLIZAMBVLFLBN"/>	valid
Confirm new password	<input type="text" value="MKLIZAMBVLXYZA"/>	unmatch ←

- Click the **Update** button in the lower right corner. Wait while the passcode is updated.



- When passcode updating is complete, the digitizer will be automatically reset. When resetting is complete, Passcode Editor will automatically close.






# Troubleshooting


Digitizer Diagnostics, Firmware Download and Passcode Editor have built-in problem detection and troubleshooting help. After launching, each program first determines whether it can communicate with the digitizer. If communication is not possible, the program provides a checklist to help you identify and correct problems. As you complete each step, click its checkbox.

After each step, the program again attempts to communicate with the digitizer. If communication is established, you will be able to start normal program operations. If communication is not established, the next step in the troubleshooting process will be displayed.


1. Make sure the digitizer is connected to the computer with the appropriate USB or SCSI cable.




2. Turn off the digitizer, then turn it on.




3. Wait three seconds.



4. Scan for hardware changes:
  - a. Right-click **My Computer**, then click **Manage** in the pop-up menu.
  - b. In the **Computer Management** window's left tree, click **Device Manager**.
  - c. In the right tree, click the computer's name (at the top of the tree).
  - d. In the toolbar, click **Scan for hardware changes** button.



4. In the **Computer Management** window, look for **Vidar Film Digitizer**, then click the corresponding checkbox in Digitizer Diagnostics.
5. To refresh this screen click on the icon the **Red Arrow** shown below.
6. You will find the digitizer listed either under Imaging Devices or Jungo depending on your operating system.



7. If you reach this point and do not find the Digitizer, no further troubleshooting suggestions are available. You can:
  - Cancel the operation and Restart both the computer and digitizer (recommended).
  - Click the **Contact form** button in the lower right corner to send a report to Vidar Technical Support. Continue at step 7 in the “Digitizer Diagnostics” chapter of this manual. See Screen Shots below

Digitizer Diagnostics 2.0

Model:

Serial number:

✔ System check

Requirements check

Run digitizer tasks

Done

✖ Checking film digitizer's connection, please wait . . .

○

Unable to find a VIDAR digitizer.

Exhausted digitizer connection troubleshooting options.

Troubleshooting Steps

☒ Check the digitizer USB or SCSI cable.

☒ Power the digitizer off, wait 3 seconds and then on.

☒ Wait for one minute.

☒ Try scanning for hardware changes from the Device Manager.

Do you see a VIDAR Film Digitizer ?

☐ YES

☒ NO    Try options 1 or 2 below.

1. Cycle power on the computer and digitizer and try again (recommended).

2. Click contact form button above to fill out contact form.

Send the report to MedTech for assistance.

VIDAR

systems corporation

The Image of Reliability

?

□

✕

CONTACT FORM

System check

Requirements check

Run digitizer tasks

Done

Model:

Serial number:

?

X

Contact Information

Please fill out the form below and click the Next button below. Fields in bold are mandatory.

**Digitizer Serial No.**

340016

**Contact Name**

John Customer

Company Name

Address

**Phone**

703-471-7070

(Include country and area codes)

**Email Address**

john.customer@hospital.com

Comments

I cannot get my digitizer to communicate with my host PC.  
Please help me  
John Customer|

V

IDAR

systems corporation

The Image of Reliability

<< BACK

NEXT >>

Adobe Acrobat Professional - [340016 - 1.pdf]

File Edit View Document Comments Tools Advanced Window Help

Create PDF Comment & Markup Send for Review Secure Sign Forms

Select 150% Help

Digitizer SN: 340016

## Film Digitizer Report - 3/1/2011 3:34:24 PM

NOTE: Please forward this report to VIDAR's Technical Support Group using one of the following options.

1. Email the report to [medtech@vidar.com](mailto:medtech@vidar.com)
2. Fax the report to the attention of Medtech at 703-471-7665.

Attn: Medical Technical Support  
 365 Herndon Parkway  
 Herndon, VA 20170  
 +1.703.471.7070 (Phone)  
 +1.703.471.7665 (Fax)  
 MedTech@vidar.com

### Digitizer Diagnostics 2.0

Model:  
Serial number:

System check  
Requirements check  
Run digitizer tests  
Done

checking film digitizer's connection, please wait ...

Unable to find a VIDAR digitizer.

Exhausted digitizer connection troubleshooting options.

#### Troubleshooting Steps

☒ Check the digitizer (USB or SCSI) cable.  
☒ Power the digitizer off, wait 3 seconds and then on.  
☒ Wait for one minute.  
☒ Try spending for hardware changes from the Device Manager.

Do you see a VIDAR Film Digitizer?

☐ YES  
☒ NO Try options 1 or 2 below.

1. Cycle power on the computer and digitizer and try again (recommended).  
 2. Click contact form button above to fill out contact form.  
 Send the report to MedTech for assistance.

CONTACT FORM

### Film Digitizer Information

Hardware driver config  
 Driver version 5, 5, 2, 0  
 Toolkit version  
 Model name  
 Password  
 Limited use device  
 Maximum bit depth  
 Maximum resolution  
 Hardware version  
 Kernel  
 Firmware version  
 DSP version  
 MCB version  
 Lamp type  
 Lamp module  
 Ballast time  
 Feeder sn  
 Feeder version  
 Camera board

### Film Digitizer USB Information

1 of 2

# Appendix

---

This appendix reproduces a typical report generated by Digitizer Diagnostics.

Reports are initially created in the *My Documents\Film Digitizer\Report* folder. You can move or copy them to any other location. Note that if an email application is detected on the host pc the program will attempt to generate an email automatically for the report to be emailed to VIDAR. You do not have to use this feature and can cancel it if you desire.

Report names have the following syntax:

serial number – the number of times the report was run .pdf

For example, a report created for digitizer serial number 350282 would have this name:

350282 - 1.pdf

If additional reports are run the – number will change:

For example

350282 – 2.pdf

350282 – 3.pdf etc.

See example report below on the next few pages.


Digitizer SN: 351727

## **Film Digitizer Report - 3/2/2011 8:38:39 AM**

NOTE: Please forward this report to VIDAR's Technical Support Group using one of the following options.

1. Email the report to [medtech@vidar.com](mailto:medtech@vidar.com)
2. Fax the report to the attention of Medtech at 703-471-7665.

Attn: Medical Technical Support  
 365 Herndon Parkway  
 Herndon, VA 20170  
 +1.703.471.7070 (Phone)  
 +1.703.471.7665 (Fax)  
 MedTech@vidar.com



### **Film Digitizer Information**

Hardware driver config	USB STI
Driver version	5, 5, 2, 0
Toolkit version	5.1.51.0
Model name	DiagnosticPRO Advantage (USA model)
Password	NGZOPBJOVLFLNN
Limited use device	False
Maximum bit depth	12
Maximum resolution	570
Hardware version	112
Kernel	9300
Firmware version	49.7
DSP version	6005
MCB version	0303
Lamp type	LED
Lamp module	White LED, full width, no rollers
Ballast time	Month: 2, Day: 8, Hour: 1
Feeder sn	193451
Feeder version	0415
Camera board	N1

### **Film Digitizer USB Information**

Device description	Not found
Device bus speed	N/A

### **Contact Information**

Contact name	John Customer
Company name	
Address	
Phone	703-471-7070
Email	<a href="mailto:john.customer@hospital.com">john.customer@hospital.com</a>

### **Film Digitizer Messages**

Main task	Digitizer Diagnostics
Main message	The operation completed successfully.
Status	Completed

### **Operating System**

OS name	Windows XP Service Pack 3 (Build 2600) 32-bit
---------	--

### **COMMENTS:**

### **SUMMARY:**

DIGITIZER INFO :Digitizer Info Passed  
 VERSION INFO :Version Information Passed  
 LIGHT OFF : 86 65 41  
 LIGHT ON :65535 65535 65535  
 FULL TEST :Full Test Passed  
 NORMALIZE :Normalize Passed





Scaler Passed  
 DSP Communication Passed  
 DSP Memory Passed  
 DSP Tests Passed  
 Full Test Passed

#### Light On

0000065415 - Ballast Level at Check  
 0000065415 - Ballast Level at Check  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535  
 65535 65535 65535

#### Light Off

87 65 45  
 86 65 45  
 85 65 46  
 88 65 45  
 96 65 45  
 86 65 48  
 87 65 45  
 85 65 39  
 87 65 40  
 86 65 41

#### Normalize

Channel 1 Zero Check 0  
 Channel 2 Zero Check 0  
 Channel 3 Zero Check 0  
 Stabilize Array Offset (495) Avg: 0  
 Stabilize Array Offset (479) Avg: 0  
 Stabilize Array Offset (463) Avg: 0  
 Stabilize Array Offset (447) Avg: 0  
 Stabilize Array Offset (431) Avg: 0  
 Stabilize Array Offset (415) Avg: 0  
 Stabilize Array Offset (399) Avg: 0  
 Stabilize Array Offset (383) Avg: 0  
 Stabilize Array Offset (367) Avg: 0  
 Stabilize Array Offset (351) Avg: 0  
 Stabilize Array Offset (335) Avg: 0  
 Stabilize Array Offset (319) Avg: 0  
 Stabilize Array Offset (303) Avg: 0  
 Stabilize Array Offset (287) Avg: 0  
 Stabilize Array Offset (271) Avg: 0  
 Stabilize Array Offset (255) Avg: 0  
 Stabilize Array Offset (239) Avg: 320  
 Stabilize Array Offset (223) Avg: 718  
 Stabilize Array Offset (223) Avg: 1116  
 Stabilize Array Offset (223) Avg: 1115  
 Array Offset Ch 1 (255) Min: 0, Max: 0  
 Array Offset Ch 1 (495) Min: 283, Max: 382  
 Array Offset Ch 1 (479) Min: 0, Max: 0  
 Array Offset Ch 1 (463) Min: 0, Max: 0  
 Array Offset Ch 1 (447) Min: 0, Max: 0  
 Array Offset Ch 1 (431) Min: 0, Max: 0  
 Array Offset Ch 1 (415) Min: 0, Max: 0  
 Array Offset Ch 1 (399) Min: 0, Max: 0  
 Array Offset Ch 1 (383) Min: 0, Max: 0  
 Array Offset Ch 1 (367) Min: 0, Max: 0  
 Array Offset Ch 1 (351) Min: 0, Max: 0  
 Array Offset Ch 1 (335) Min: 0, Max: 0  
 Array Offset Ch 1 (319) Min: 0, Max: 0  
 Array Offset Ch 1 (303) Min: 0, Max: 0  
 Array Offset Ch 1 (287) Min: 0, Max: 0  
 Array Offset Ch 1 (271) Min: 0, Max: 0  
 Array Offset Ch 1 (270) Min: 0, Max: 25  
 Array Offset Ch 1 (269) Min: 0, Max: 47  
 Array Offset Ch 1 (268) Min: 0, Max: 66  
 Array Offset Ch 1 (267) Min: 0, Max: 88  
 Array Offset Ch 1 (266) Min: 21, Max: 114  
 Array Offset Ch 1 (265) Min: 40, Max: 137  
 Array Offset Ch 1 (264) Min: 67, Max: 172  
 Array Offset Ch 1 (263) Min: 88, Max: 184  
 Array Offset Ch 1 (262) Min: 121, Max: 221  
 Array Offset Ch 1 (262) Min: 145, Max: 239

Array Offset Ch 2 (255) Min: 0, Max: 0  
 Array Offset Ch 2 (495) Min: 147, Max: 240  
 Array Offset Ch 2 (479) Min: 0, Max: 0  
 Array Offset Ch 2 (463) Min: 0, Max: 0  
 Array Offset Ch 2 (447) Min: 0, Max: 0  
 Array Offset Ch 2 (431) Min: 0, Max: 0  
 Array Offset Ch 2 (415) Min: 0, Max: 0  
 Array Offset Ch 2 (399) Min: 0, Max: 0  
 Array Offset Ch 2 (383) Min: 0, Max: 0  
 Array Offset Ch 2 (367) Min: 0, Max: 0  
 Array Offset Ch 2 (351) Min: 0, Max: 0  
 Array Offset Ch 2 (335) Min: 0, Max: 0  
 Array Offset Ch 2 (319) Min: 0, Max: 0  
 Array Offset Ch 2 (303) Min: 0, Max: 0  
 Array Offset Ch 2 (287) Min: 0, Max: 0  
 Array Offset Ch 2 (271) Min: 0, Max: 0  
 Array Offset Ch 2 (270) Min: 0, Max: 0  
 Array Offset Ch 2 (269) Min: 0, Max: 0  
 Array Offset Ch 2 (268) Min: 0, Max: 0  
 Array Offset Ch 2 (267) Min: 0, Max: 0  
 Array Offset Ch 2 (266) Min: 0, Max: 0  
 Array Offset Ch 2 (265) Min: 0, Max: 6  
 Array Offset Ch 2 (264) Min: 0, Max: 23  
 Array Offset Ch 2 (263) Min: 0, Max: 47  
 Array Offset Ch 2 (262) Min: 0, Max: 76  
 Array Offset Ch 2 (261) Min: 4, Max: 98  
 Array Offset Ch 2 (260) Min: 28, Max: 124  
 Array Offset Ch 2 (259) Min: 49, Max: 151  
 Array Offset Ch 2 (258) Min: 82, Max: 175  
 Array Offset Ch 2 (257) Min: 101, Max: 205  
 Array Offset Ch 2 (256) Min: 128, Max: 227  
 Array Offset Ch 2 (256) Min: 144, Max: 237  
 Array Offset Ch 3 (255) Min: 0, Max: 0  
 Array Offset Ch 3 (495) Min: 174, Max: 255  
 Array Offset Ch 3 (479) Min: 0, Max: 0  
 Array Offset Ch 3 (463) Min: 0, Max: 0  
 Array Offset Ch 3 (447) Min: 0, Max: 0  
 Array Offset Ch 3 (431) Min: 0, Max: 0  
 Array Offset Ch 3 (415) Min: 0, Max: 0  
 Array Offset Ch 3 (399) Min: 0, Max: 0  
 Array Offset Ch 3 (383) Min: 0, Max: 0  
 Array Offset Ch 3 (367) Min: 0, Max: 0  
 Array Offset Ch 3 (351) Min: 0, Max: 0  
 Array Offset Ch 3 (335) Min: 0, Max: 0  
 Array Offset Ch 3 (319) Min: 0, Max: 0  
 Array Offset Ch 3 (303) Min: 0, Max: 0  
 Array Offset Ch 3 (287) Min: 0, Max: 0  
 Array Offset Ch 3 (271) Min: 0, Max: 0  
 Array Offset Ch 3 (270) Min: 0, Max: 0  
 Array Offset Ch 3 (269) Min: 0, Max: 0  
 Array Offset Ch 3 (268) Min: 0, Max: 0  
 Array Offset Ch 3 (267) Min: 0, Max: 0  
 Array Offset Ch 3 (266) Min: 0, Max: 0  
 Array Offset Ch 3 (265) Min: 0, Max: 16  
 Array Offset Ch 3 (264) Min: 0, Max: 55  
 Array Offset Ch 3 (263) Min: 0, Max: 63  
 Array Offset Ch 3 (262) Min: 11, Max: 102  
 Array Offset Ch 3 (261) Min: 32, Max: 116  
 Array Offset Ch 3 (260) Min: 55, Max: 145  
 Array Offset Ch 3 (259) Min: 78, Max: 162  
 Array Offset Ch 3 (258) Min: 108, Max: 193  
 Array Offset Ch 3 (257) Min: 126, Max: 212  
 Array Offset Ch 3 (257) Min: 153, Max: 237  
 Pixel Offset (Ch 1) (Target 1) Min: 174, Avg: 179, Max: 223  
 Pixel Offset (Ch 2) (Target 1) Min: 178, Avg: 182, Max: 225  
 Pixel Offset (Ch 3) (Target 1) Min: 184, Avg: 188, Max: 222  
 Elapsed Time 62.730  
 Pixel Offset Check (Ch 1) Min: 2, Avg: 3, Max: 5  
 Pixel Offset Check (Ch 2) Min: 2, Avg: 3, Max: 5  
 Pixel Offset Check (Ch 3) Min: 2, Avg: 3, Max: 5  
 Elapsed Time 67.350  
 Channel 1 Ones Check 65535  
 Channel 2 Ones Check 65535  
 Channel 3 Ones Check 65535  
 Ballast Find High Channel 1 65362  
 Ballast Find High Channel 2 65358  
 Ballast Find High Channel 3 65352  
 Channel Selected 1  
 Elapsed Time 72.590  
 Ballast Adjust (230) Min: 61372, Max: 65362  
 Video LUT Factor 00001.0000  
 Elapsed Time 72.845  
 Ballast Adjust (229) Ch 1 Min: 61384, Max: 65362  
 Ballast Adjust (228) Ch 1 Min: 61114, Max: 65362  
 Ballast Adjust (227) Ch 1 Min: 60701, Max: 65362

Digitizer SN: 351727

Ballast Adjust (226) Ch 1 Min: 60648, Max: 65362  
 Ballast Adjust (225) Ch 1 Min: 60314, Max: 65362  
 Ballast Adjust (224) Ch 1 Min: 60070, Max: 65362  
 Ballast Adjust (223) Ch 1 Min: 59895, Max: 65362  
 Ballast Adjust (222) Ch 1 Min: 59735, Max: 65362  
 Ballast Adjust (221) Ch 1 Min: 59382, Max: 65362  
 Ballast Adjust (220) Ch 1 Min: 59025, Max: 65362  
 Ballast Adjust (219) Ch 1 Min: 58900, Max: 65362  
 Ballast Adjust (218) Ch 1 Min: 58605, Max: 65362  
 Ballast Adjust (217) Ch 1 Min: 58327, Max: 65362  
 Ballast Adjust (216) Ch 1 Min: 58023, Max: 65362  
 Ballast Adjust (215) Ch 1 Min: 57861, Max: 65362  
 Ballast Adjust (214) Ch 1 Min: 57592, Max: 65362  
 Ballast Adjust (213) Ch 1 Min: 57342, Max: 65362  
 Ballast Adjust (212) Ch 1 Min: 57008, Max: 65362  
 Ballast Adjust (211) Ch 1 Min: 56729, Max: 65362  
 Ballast Adjust (210) Ch 1 Min: 56567, Max: 65362  
 Ballast Adjust (209) Ch 1 Min: 56271, Max: 65362  
 Ballast Adjust (208) Ch 1 Min: 55994, Max: 65362  
 Ballast Adjust (207) Ch 1 Min: 55686, Max: 65362  
 Ballast Adjust (206) Ch 1 Min: 55586, Max: 65362  
 Ballast Adjust (205) Ch 1 Min: 55151, Max: 65362  
 Ballast Adjust (204) Ch 1 Min: 54831, Max: 65362  
 Ballast Adjust (203) Ch 1 Min: 54717, Max: 65362  
 Ballast Adjust (202) Ch 1 Min: 54529, Max: 65362  
 Ballast Adjust (201) Ch 1 Min: 54239, Max: 65362  
 Ballast Adjust (200) Ch 1 Min: 53839, Max: 65362  
 Ballast Adjust (199) Ch 1 Min: 53812, Max: 65362  
 Ballast Adjust (198) Ch 1 Min: 53349, Max: 65362  
 Ballast Adjust (197) Ch 1 Min: 53083, Max: 65362  
 Ballast Adjust (196) Ch 1 Min: 52859, Max: 65362  
 Ballast Adjust (195) Ch 1 Min: 52588, Max: 65362  
 Ballast Adjust (194) Ch 1 Min: 52364, Max: 65362  
 Ballast Adjust (193) Ch 1 Min: 52154, Max: 65362  
 Ballast Adjust (192) Ch 1 Min: 51728, Max: 65362  
 Ballast Adjust (191) Ch 1 Min: 51561, Max: 65362  
 Ballast Adjust (190) Ch 1 Min: 51381, Max: 65362  
 Ballast Adjust (189) Ch 1 Min: 51134, Max: 65362  
 Ballast Adjust (188) Ch 1 Min: 50887, Max: 65362  
 Ballast Adjust (187) Ch 1 Min: 50457, Max: 65362  
 Ballast Adjust (186) Ch 1 Min: 50255, Max: 65362  
 Ballast Adjust (185) Ch 1 Min: 49908, Max: 65362  
 Ballast Adjust (184) Ch 1 Min: 49795, Max: 65362  
 Ballast Adjust (183) Ch 1 Min: 49547, Max: 65362  
 Ballast Adjust (182) Ch 1 Min: 49287, Max: 65362  
 Ballast Adjust (181) Ch 1 Min: 49044, Max: 65362  
 Ballast Adjust (180) Ch 1 Min: 48672, Max: 65362  
 Ballast Adjust (179) Ch 1 Min: 48529, Max: 65362  
 Ballast Adjust (178) Ch 1 Min: 48221, Max: 65362  
 Ballast Adjust (177) Ch 1 Min: 47966, Max: 65362  
 Ballast Adjust (176) Ch 1 Min: 47612, Max: 65362  
 Ballast Adjust (175) Ch 1 Min: 47329, Max: 65362  
 Ballast Adjust (174) Ch 1 Min: 47058, Max: 65362  
 Ballast Adjust (173) Ch 1 Min: 46881, Max: 65362  
 Ballast Adjust (172) Ch 1 Min: 46586, Max: 65362  
 Ballast Adjust (171) Ch 1 Min: 46381, Max: 65362  
 Ballast Adjust (170) Ch 1 Min: 46035, Max: 65362  
 Ballast Adjust (169) Ch 1 Min: 45912, Max: 65362  
 Ballast Adjust (168) Ch 1 Min: 45591, Max: 65362  
 Ballast Adjust (167) Ch 1 Min: 45399, Max: 65362  
 Ballast Adjust (166) Ch 1 Min: 44994, Max: 65362  
 Ballast Adjust (165) Ch 1 Min: 44704, Max: 65362  
 Ballast Adjust (164) Ch 1 Min: 44534, Max: 65362  
 Ballast Adjust (163) Ch 1 Min: 44223, Max: 65362  
 Ballast Adjust (162) Ch 1 Min: 44057, Max: 65362  
 Ballast Adjust (161) Ch 1 Min: 43744, Max: 65362  
 Ballast Adjust (160) Ch 1 Min: 43433, Max: 65362  
 Ballast Adjust (159) Ch 1 Min: 43285, Max: 65362  
 Ballast Adjust (158) Ch 1 Min: 42915, Max: 65362  
 Ballast Adjust (157) Ch 1 Min: 42663, Max: 65362  
 Ballast Adjust (156) Ch 1 Min: 42379, Max: 65361  
 Ballast Adjust (155) Ch 1 Min: 42114, Max: 65361  
 Ballast Adjust (154) Ch 1 Min: 41908, Max: 65361  
 Ballast Adjust (153) Ch 1 Min: 41655, Max: 65361  
 Ballast Adjust (152) Ch 1 Min: 41337, Max: 65361  
 Ballast Adjust (151) Ch 1 Min: 40986, Max: 65361  
 Ballast Adjust (150) Ch 1 Min: 40753, Max: 65360  
 Ballast Adjust (149) Ch 1 Min: 40550, Max: 65359  
 Ballast Adjust (148) Ch 1 Min: 40222, Max: 64921  
 Ballast Adjust (147) Ch 1 Min: 39893, Max: 64521  
 Ballast Adjust (146) Ch 1 Min: 39646, Max: 64113  
 Ballast Adjust (145) Ch 1 Min: 39425, Max: 63725  
 Ballast Adjust (144) Ch 1 Min: 39091, Max: 63238  
 Ballast Adjust (143) Ch 1 Min: 38884, Max: 62909  
 Ballast Adjust (142) Ch 1 Min: 38630, Max: 62485

Ballast Adjust (141) Ch 1 Min: 38333, Max: 62080  
 Ballast Adjust (140) Ch 1 Min: 38092, Max: 61760  
 Ballast Adjust (140) Ch 1 Min: 37861, Max: 61304  
 Ballast Adjust (140) Ch 2 Min: 42450, Max: 61154  
 Ballast Adjust (140) Ch 3 Min: 37962, Max: 59547  
 Elapsed Time 96.485  
 Pixel limited Left 0 Right 1 on Ch 1  
 Pixel Gain (Ch 1) Min (Gain): 37964 (41492), Max (Gain): 61091 ( 974)  
 Pixel Gain (Limited) Min (Gain): 37964 (41492), Max (Gain): 61091 ( 974)  
 Pixel Gain (Ch 1) Low Count (Min): 0 (61959), High Count (Max): 0 (62037) of 7780  
 Pixel limited Left 0 Right 1 on Ch 2  
 Pixel Gain (Ch 2) Min (Gain): 42620 (29800), Max (Gain): 60974 ( 1102)  
 Pixel Gain (Limited) Min (Gain): 42961 (29042), Max (Gain): 60974 ( 1102)  
 Pixel Gain (Ch 2) Low Count (Min): 0 (61953), High Count (Max): 0 (62040) of 7780  
 Pixel limited Left 0 Right 1 on Ch 3  
 Pixel Gain (Ch 3) Min (Gain): 38064 (41211), Max (Gain): 59377 ( 2895)  
 Pixel Gain (Limited) Min (Gain): 38064 (41211), Max (Gain): 59377 ( 2895)  
 Pixel Gain (Ch 3) Low Count (Min): 0 (61957), High Count (Max): 0 (62042) of 7780  
 Elapsed Time 105.235  
 Second Offset Correction  
 Elapsed Time 105.430  
 Summary Information  
 Channel Offset 00262 00256 00257  
 Pixel Offset Min 00174 00178 00184  
 Pixel Offset Ave 00179 00182 00188  
 Pixel Offset Max 00223 00225 00222  
 Ballast Level 140  
 Video LUT Factor 00001.0000  
 Pixel Gain Min 61959 61953 61957  
 Pixel Gain Max 62037 62040 62042  
 Elapsed Time 105.445  
 Normalize Passed

Digitizer SN: 351727